

iPS cell induction from human non-T, B cells from peripheral blood

Keisuke Okita

<day 0>

1. Culture medium preparation:

- a. Add 2 mL of the following culture medium containing cytokines to a well 6 well plate.

StemSpan H3000
100 ng/mL IL-6
300 ng/mL SCF
300 ng/mL TPO
300 ng/mL Flt3 ligand
10 ng/mL IL-3

- b. Add 2 mL of PBS in the remaining wells of the plate and store the plate at 37°C, 5% CO₂.

2. Use Ficoll-paque Plus to purify mononuclear cells:

- a. To 10 mL of anti-coagulated blood (EDTA), add 10 ml of PBS.
- b. In two 15 mL tubes, add 5 mL Ficoll-Paque and the gently add 10 mL of blood+PBS
- c. Spin tubes at 400 X g for 30 mins. at 18°C. Use slow acceleration and no brake.
- d. Remove the plasma for the top fraction without disrupting the mononuclear cells at the interface.
- e. Transfer the cells at interface to a 15 mL tube and add 12 mL of PBS.
- f. Spin at 200 X g for 10 mins. at 18°C (no brake).
- g. Resuspend cells in 3 mL of H3000 medium and count cells.
- h. Prepare aliquots of 3 x 10⁶ cells in 1.5 mL tube.
- i. Spin at 200 X g for 10 mins. at 18°C (no brake).
- j. Aspirate supernatant.

3. Plating:

- a. Resuspend the cells in the medium prepared in 1.
- b. Store the plate at 37°C, 5% CO₂ for around 6 days. Medium change is not needed.

<day 5>

1. Prepare 6 well plate coated with MEF feeder cells. (3 X 10⁵ cell/well)
2. iPSCs can be established in non-feeder condition, but the efficiency is low.
 - a. Add 2 ml of RetroNectin solution (20 µg/ml) to a well 6 well plate .
 - b. Store at RT for 2 hr.
 - c. Wash the plate once with PBS before cell seeding.

<day 6>

1. Culture medium preparation:

- a. Prepare the following culture medium.

StemSpan H3000
100 ng/mL IL-6
300 ng/mL SCF
300 ng/mL TPO
300 ng/mL Flt3 ligand
10 ng/mL IL-3

2. Harvest the cultured cells:

- a. Resuspend the cells in the medium and count cells. Number of live cells is usually around 1×10^6 .
- b. Harvest the floating cells into 15 ml tube.
- c. Spin at 200 X g for 10 mins. at 18°C (no brake). During spin, prepare nucleofection solution.

3. Nucleofection:

Amaya CD34 Solution	81.8 μ L
Supplement	18.2 μ L
Plasmid mix	3 μ L (3 μ g)

- a. Aspirate supernatant completely by hand using a pipette.
- b. Add nucleofection solution and suspend cells, be careful not to create any bubbles.
- c. Perform nucleofection using pre-stored program U-008.

4. Plating, 10^6 to 10^4 cells per well:

- a. Immediately following nucleofection, add 800 μ L of H3000 to the electroporation cuvette, and harvest the cells. Metal ions in the nucleofection solution are toxic to cells!
- b. Plate the cells to MEF feeder or RetroNectin coated 6-well plate ranging from 5×10^5 cells to 1×10^4 cells per well with the medium prepared above.

<day 8, 10, and 12>

Add additional 1.5 mL of TeSR2 medium per well.

<day 14 - >

Replace medium with 1.5 mL of TeSR2 per well.

* Medium replacement is performed every 2 days.

<day 25 to 35>

Pick colonies of about 2 mm diameter.

Reagents

Ficoll-paque Plus (GE Healthcare, 17-1440-02)

StemSpan H3000 (StemCell Technologies, 09800)

IL-6 (PeproTech, AF-200-06B)

SCF (PeproTech, AF-300-07B)

TPO (PeproTech, AF-300-18B)

Flt3 ligand (PeproTech, AF-300-19B)
 IL-3 (PeproTech, AF-200-03B)
 Amaxa Human CD34+ cell Nucleofector Kit (Lonza, VPA-1003)
 Nucleofector 2b (Lonza, AAB-1001)
 RetroNectin (Takara, T100A)
 TeSR2 (StemCell Technologies, ST-05860)
 Plasmid (Addgene, http://www.addgene.org/Shinya_Yamanaka)

Use following plasmid mixtures. Set 1 shows high efficiency. In set 2, we omitted WPRE sequence and replaced shRNA against p53 with dominant negative form of mouse p53, which exist in set 1.

Plasmid set 1	pCXLE-hOCT3/4-shp53-F	0.83 µg
	pCXLE-hSK	0.83 µg
	pCXLE-hUL	0.83 µg
	pCXWB-EBNA1	0.5 µg
Plasmid set 2	pCE-hOCT3/4	0.63 µg
	pCE-hSK	0.63 µg
	pCE-hUL	0.63 µg
	pCE-mp53DD	0.63 µg
	pCXB-EBNA1	0.5 µg

Reference

1. Okita, K., et al. An Efficient Non-viral Method to Generate Integration-Free Human iPS Cells from Cord Blood and Peripheral Blood Cells. *Stem Cells*. 2012 Nov 29
2. Okita, K., et al. A more efficient method to generate integration-free human iPS cells. *Nat Methods*. 8, 409-12 (2011).
3. Mack, AA., et al. Generation of Induced Pluripotent Stem Cells from CD34+ Cells across Blood Drawn from Multiple Donors with Non-Integrating Episomal Vectors. *PLoS One*. 6, e27956 (2011).